



## Human CD218a-IL18R1 recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P6206
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Met1-Arg329
<b>TAG</b>	C-Terminal Fc Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	64.23 kDa
<b>FORM</b>	
<b>BUFFER</b>	0.01M PBS, pH 7.4.
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	Q13478
<b>ALIAS/SYNONYMS</b>	CD218a, IL-18R1, IL18R1, CDw218a, IL-18R-alpha, IL-18R-1, Interleukin-18 receptor 1, IL-18Ralpha, IL1 receptor-related protein, CD218 antigen-like family member A, IL-1Rrp, IL1RRP, IL1R-rp, Interleukin-18 receptor alpha

