



Human CD111-NECTIN1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6207
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Thr334
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	38.02 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q15223
ALIAS/SYNONYMS	Herpesvirus entry mediator C, CD111, Herpes virus entry mediator C, HIgR, Poliovirus receptor-related protein 1, PRR1, NECTIN1, Nectin cell adhesion molecule 1, Nectin-1, HveC, HVEC, Herpesvirus Ig-like receptor, PVRL1

