



Human AMHR2 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6208
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Ser144
TAG	C-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	41.65 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q16671
ALIAS/SYNONYMS	MIS type II receptor, MISR2, AMHR2, AMHR, Anti-Muellerian hormone type II receptor, MRII, AMH type II receptor, MISRII, Anti-Muellerian hormone type-2 receptor

