



Human AOC3-VAP-1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6209
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Gly27-Asn763
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Insect
MOLECULAR WEIGHT	84.75 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q16853
ALIAS/SYNONYMS	Membrane primary amine oxidase, AOC3, SSAO, Semicarbazide-sensitive amine oxidase, Vascular adhesion protein 1, VAP-1, HPAO, Copper amine oxidase, VAP1

