



Human CD300g-CD300LG recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6213
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Leu249
TAG	C-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	53.41 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q6UXG3
ALIAS/SYNONYMS	CD300g, CLM-9, Triggering receptor expressed on myeloid cells 4, TREM-4, TREM4, CLM9, CD300LG, CD300 antigen-like family member G, CMRF35-like molecule 9

