



Human CD109 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6214
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Val22-Ser1268
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	142.98 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q6YHK3
ALIAS/SYNONYMS	CPAMD7, r150, CD109, 150 kDa TGF-beta-1-binding protein, Platelet-specific Gv antigen, CD109 antigen, p180, C3 and PZP-like alpha-2-macroglobulin domain-containing protein 7

