



Human CD85h-LILRA2 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6216
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Asn449
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	50.34 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8N149
ALIAS/SYNONYMS	Leukocyte immunoglobulin-like receptor subfamily A member 2, LIR-7, LILRA2, LIR7, Leukocyte immunoglobulin-like receptor 7, CD85h, Immunoglobulin-like transcript 1, ILT1, CD85 antigen-like family member H, ILT-1

