



Human CD85d-LILRB2/ILT4 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6217
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Val461
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	50.88 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8N423
ALIAS/SYNONYMS	CD85 antigen-like family member D, Monocyte/macrophage immunoglobulin-like receptor 10, MIR-10, ILT4, LIR-2, ILT-4, Immunoglobulin-like transcript 4, MIR10, CD85d, Leukocyte immunoglobulin-like receptor subfamily B member 2, Leukocyte immunoglobulin-like

