



Human CD85e-LILRA3 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6218
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Glu439
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	48.46 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8N6C8
ALIAS/SYNONYMS	Immunoglobulin-like transcript 6, ILT-6, CD85e, LILRA3, Monocyte inhibitory receptor HM43/HM31, LIR-4, CD85 antigen-like family member E, ILT6, Leukocyte immunoglobulin-like receptor 4, LIR4, Leukocyte immunoglobulin-like receptor subfamily A member 3

