



Human CD85j-LILRB1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6220
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-His458
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	50.59 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8NHL6
ALIAS/SYNONYMS	LIR1, Immunoglobulin-like transcript 2, MIR7, CD85 antigen-like family member J, Monocyte/macrophage immunoglobulin-like receptor 7, ILT-2, Leukocyte immunoglobulin-like receptor subfamily B member 1, LIR-1, CD85j, ILT2, LILRB1, MIR-7, Leukocyte immunoglo

