



Human CD352-SLAMF6 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6223
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Met226
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	26.43 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q96DU3
ALIAS/SYNONYMS	CD352, NTB-A, NK-T-B-antigen, KALI, SLAMF6, Activating NK receptor, SLAM family member 6

