



Human LINGO1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6224
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Cys42-Pro549
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	60.33 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q96FE5
ALIAS/SYNONYMS	LRRN6A, LERN1, Leucine-rich repeat neuronal protein 1, Leucine-rich repeat and immunoglobulin domain-containing protein 1, LINGO1, Leucine-rich repeat neuronal protein 6A, Leucine-rich repeat and immunoglobulin-like domain-containing nogo receptor-interac

