



## Human CD307a-FCRL1 recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P6225
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Met1-Asn303
<b>TAG</b>	C-Terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	33.94 kDa
<b>FORM</b>	
<b>BUFFER</b>	0.01M PBS, pH 7.4.
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	Q96LA6
<b>ALIAS/SYNONYMS</b>	CD307a, Fc receptor homolog 1, hIFGP1, FcRL1, Fc receptor-like protein 1, FCRL1, IRTA5, IFGP family protein 1, FcR-like protein 1, FcRH1, Immune receptor translocation-associated protein 5, IFGP1, FCRH1

