



Human CD354-TREM1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6227
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Arg200
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	23.72 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9NP99
ALIAS/SYNONYMS	TREM1, Triggering receptor expressed on monocytes 1, Triggering receptor expressed on myeloid cells 1, TREM-1, CD354

