



Human CD320 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6228
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Gly208
TAG	C-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	47.87 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9NPF0
ALIAS/SYNONYMS	FDC-SM-8D6, FDC-signaling molecule 8D6, CD320, TCb1R, 8D6 antigen, 8D6A, Transcobalamin receptor, CD320 antigen

