



Human CD207 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6232
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Pro65-Pro328
TAG	N-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	57.40 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9UJ71
ALIAS/SYNONYMS	C-type lectin domain family 4 member K, CLEC4K, Langerin, CD207

