



Human CD201-PROCR recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6233
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Ser210
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	24.86 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9UNN8
ALIAS/SYNONYMS	Endothelial protein C receptor, PROCR, EPCR, APC receptor, Endothelial cell protein C receptor, CD201, Activated protein C receptor

