



Human CD208-LAMP3 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6234
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Thr381
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	41.52 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9UQV4
ALIAS/SYNONYMS	DCLAMP, CD208, LAMP-3, LAMP3, Lysosomal-associated membrane protein 3, TSC403, Lysosome-associated membrane glycoprotein 3, DC LAMP, Protein TSC403, DC-lysosome-associated membrane glycoprotein

