



Human CD307e-FCRL5 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6235
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Arg844
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	93.24 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q96RD9
ALIAS/SYNONYMS	FcRL5, FCRH5, FCRL5, Fc receptor homolog 5, FcRH5, CD307e, FcR-like protein 5, IRTA2, Immune receptor translocation-associated protein 2, BXMAS1, Fc receptor-like protein 5

