



## Human CD367-CLEC4A recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P6237
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Gln70-Leu237
<b>TAG</b>	C-Terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	23.05 kDa
<b>FORM</b>	
<b>BUFFER</b>	0.01M PBS, pH 7.4.
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	Q9UMR7
<b>ALIAS/SYNONYMS</b>	C-type lectin DDB27, CLECSF6, CLEC4A, LLIR, CD367, Lectin-like immunoreceptor, DCIR, C-type lectin domain family 4 member A, Dendritic cell immunoreceptor, C-type lectin superfamily member 6

