



Human ANGPT2-Ang2 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6245
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Tyr19-Phe496
TAG	C-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	83.24 kDa
FORM	
BUFFER	Lyophilized from a solution in PBS pH 7.4, 1mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O15123
ALIAS/SYNONYMS	Angiopoietin-2, ANG-2, ANGPT2

