



Human IGF1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6251
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Ala118
TAG	C-Terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	39.4 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P05019
ALIAS/SYNONYMS	IGF-I, IGF1, Mechano growth factor, MGF, Insulin-like growth factor I, IBP1, Somatomedin-C

