



Human OX40L recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6269
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Gln51-Leu183
TAG	N-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	14.63 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P23510
ALIAS/SYNONYMS	TXGP1, TNFSF4

