



Human IL12b and IL23a recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6270
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	(Arg20-Pro189) & (Ile23-Ser328)
TAG	N-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	Reduced: 20 kDa & 44 kDa; Nonreduced: 60 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9NPF7 & P29460
ALIAS/SYNONYMS	Interleukin-23, IL23, IL-23, interleukin 12

