



Human CD20/MS4A1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6274
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ile143-Ser188
TAG	N-Terminal His & N-Terminal Trx Tags
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	19.58 kDa
FORM	
BUFFER	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P11836
ALIAS/SYNONYMS	B-lymphocyte surface antigen B1, Membrane-spanning 4-domains subfamily A member 1, MS4A1, Leukocyte surface antigen Leu-16, B-lymphocyte antigen CD20, Bp35, CD20

