



Human IL20RB recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6279
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ala37-Glu290
TAG	N-Terminal His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	30.82 kDa
BUFFER	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	4 °C for short term (1 week), -20 °C or -80 °C for long term (avoid freezing/thawing cycles; the addition of 20 - 40% glycerol improves cryoprotection)
UNIPROT	Q6UXL0
ALIAS/SYNONYMS	IL-20R-beta, IL-20R2, FNDC6, IL-20RB, DIRS1, IL20RB, Interleukin-20 receptor subunit beta, Fibronectin type III domain containing 6, IL-20 receptor subunit beta

