



Human IL26/AK155 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6281
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Lys22-Gln171
TAG	C-Terminal His Tag

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	18.71 kDa
BUFFER	Lyophilized from a solution in PBS pH 7.4, 0.02% NLS, 1mM EDTA, 4% Trehalose, 1% Mannitol.
STABILITY & STORAGE	4 °C for short term (1 week), -20 °C or -80 °C for long term (avoid freezing/thawing cycles; the addition of 20 - 40% glycerol improves cryoprotection)
UNIPROT	Q9NPH9
ALIAS/SYNONYMS	IL-26, IL26, Protein AK155, Interleukin-26, AK155

