



Ustekinumab Biosimilar - Anti-Human IL-12 IL-23 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1011
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
SOURCE	DrugBank DB05679

PRODUCT INFORMATION

ORIGIN SPECIES	Human
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS pH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Ustekinumab,Ustekinumab,Human IL-12 IL-23,anti-Human IL-12 IL-23

