



Utomilumab Biosimilar - Anti-Human TNFRSF9 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1012
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2
SOURCE	DrugBank DB15113

PRODUCT INFORMATION

ORIGIN SPECIES	Human
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS pH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Utomilumab,Utomilumab,Human TNFRSF9,anti-Human TNFRSF9

