



Zolbetuximab Biosimilar - Anti-CLDN18.2 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1013
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
SOURCE	DrugBank DB15118

PRODUCT INFORMATION

ORIGIN SPECIES	Human
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS pH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Zolbetuximab,Zolbetuximab,CLDN18.2,anti-CLDN18.2

