



Dulaglutide Biosimilar - Anti-Glucagon-like peptide-1 receptor agonist mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1017
CLONALITY	Monoclonal Antibody
ISOTYPE	hG4-Fc
SOURCE	DrugBank DB09045

PRODUCT INFORMATION

ORIGIN SPECIES	Human
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS pH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Dulaglutide,Dulaglutide,Glucagon-like peptide-1 receptor agonist,anti-Glucagon-like peptide-1 receptor agonist

