



Certolizumab Biosimilar - Anti-TNF-Alpha mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1024
CLONALITY	Monoclonal Antibody
ISOTYPE	Fab'G1Kappa
SOURCE	DrugBank DB08904

PRODUCT INFORMATION

ORIGIN SPECIES	Human
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS pH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Certolizumab,Certolizumab,TNF-Alpha,anti-TNF-Alpha

