



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Romilkimab Biosimilar - Anti-IL13, IL4 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1027
CLONALITY	Monoclonal Antibody
ISOTYPE	[VH-VH' Gamma1 + VL-VL' kappa (bispecific)
SOURCE	IMGT/mAb-DB 794

PRODUCT INFORMATION

ORIGIN SPECIES	Human
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS pH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Romilkimab,Romilkimab,IL13, IL4,anti-IL13, IL4



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.