



Epratuzumab Biosimilar - Anti-CD22 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1042
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-nd
SOURCE	CAS 205923-57-5

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Epratuzumab,90Y-DOTA-hLL2,90Y-epratuzumab tetraxetan,AMG-412,90Y-hLL2,hLL2,CD22,anti-CD22

