



Gemtuzumab Biosimilar - Anti-CD33 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1044
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-kappa
SOURCE	CAS 220578-59-6

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Gemtuzumab,CDP-771,CMA-676 (ADC),WAY-CMA-676,gemtuzumab ozogamicin,hP67.6 (naked),hP67.6-calicheamicin (ADC),CD33,anti-CD33

