



## Blinatumomab Biosimilar - Anti-CD19, CD3E mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1067
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	(scFv-kappa-heavy)-(scFv-heavy-kappa); His-Tag
<b>SOURCE</b>	CAS 853426-35-4

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Mus musculus
<b>EXPRESSION SYSTEM</b>	Mammalian cells
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	Blinatumomab,AMG103,BITE MT-103,MEDI-538,MT103,bscCD19xCD3,CD19, CD3E,anti-CD19, CD3E

