



Abagovomab Biosimilar - Anti-idiotope of anti-Mus musculus monoclonal antibody OC126 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1075
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS 792921-10-9

PRODUCT INFORMATION

ORIGIN SPECIES	Mus musculus
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Abagovomab,ACA125,MEN-2234,idiotope of anti-Mus musculus monoclonal antibody OC126,anti-idiotope of anti-Mus musculus monoclonal antibody OC126

