



Imciromab Biosimilar - Anti-Myosin, Cardiac mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1085
CLONALITY	Monoclonal Antibody
ISOTYPE	Fab-G2a-kappa
SOURCE	CAS 126132-83-0

PRODUCT INFORMATION

ORIGIN SPECIES	Mus musculus
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Imciromab,R1D10 pentetate,Myosin, Cardiac,anti-Myosin, Cardiac

