



Biciromab Biosimilar - Anti-Fibrin II beta chain, fibrin mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1091
CLONALITY	Monoclonal Antibody
ISOTYPE	Fab'-G1-kappa
SOURCE	CAS 138783-13-8

PRODUCT INFORMATION

ORIGIN SPECIES	Mus musculus
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Biciromab, Indium (111In) biciromab, T2G1, Fibrin II beta chain, fibrin, anti-Fibrin II beta chain, fibrin

