



## 12D10-Fab Biosimilar - Anti-F7 mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1110
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	Fab-nd-nd
<b>SOURCE</b>	

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Chimeric
<b>EXPRESSION SYSTEM</b>	Mammalian cells
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	12D10-Fab,12D10 Fab,12D10-Fab,F7,anti-F7

