



Lexatumumab Biosimilar - Anti-TNFRSF10B, CD262 mAb - Research Grade

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PX-TA1135 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1-lambda |
| SOURCE | CAS 845816-02-6 |

PRODUCT INFORMATION

| | |
|---------------------------|---|
| ORIGIN SPECIES | Homo sapiens |
| EXPRESSION SYSTEM | Mammalian cells |
| BUFFER | PBS buffer PH7.5 |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |
| ALIAS/SYNONYMS | Lexatumumab,HGS-ETR2,TNFRSF10B, CD262,anti-TNFRSF10B, CD262 |

