



Ofatumumab Biosimilar - Anti-MS4A1, CD20 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1146
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS 679818-59-8

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Ofatumumab,HuMax-CD20,MS4A1, CD20,anti-MS4A1, CD20

