



# ProteoGenix

## ProteoGenix S.A.S

19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70

Fax +33 (0) 978 53 36 90

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## ProteoGenix INC.

1000 Bearcat Way, Suite  
100, Morrisville, NC 27560,  
USA

Tel +1(919)234-1277

[products@proteogenix.fr](mailto:products@proteogenix.fr)

## PRODUCT DATASHEET

### OVERVIEW

---

<b>PRODUCT NAME</b>	Panitumumab Biosimilar - Anti-EGFR mAb - Research Grade
<b>PRODUCT REFERENCE</b>	PX-TA1148
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG2, Kappa
<b>SOURCE</b>	CAS 339177-26-3

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Homo sapiens
<b>PRODUCTION HOST</b>	Mammalian cells
<b>FORM</b>	Liquid
<b>BUFFER</b>	PBS buffer PH7.5
<b>STABILITY &amp; STORAGE</b>	Store at 4°C for short term. Freeze and store at -20°C or -80°C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Avoid freeze/thaw cycles, aliquot the product according to your needs.

You can contact us at [products@proteogenix.fr](mailto:products@proteogenix.fr) for any question

**OUR PRODUCTS ARE FOR RESEARCH USE ONLY**