



Dorlimomab Biosimilar - Anti-(ricin A chain) immunotoxin mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1155
CLONALITY	Monoclonal Antibody
ISOTYPE	F(ab') ₂ -nd-nd
SOURCE	

PRODUCT INFORMATION

ORIGIN SPECIES	Mus musculus
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Dorlimomab,4197X-RA,MDX-RA (ricin A chain) immunotoxin,(ricin A chain) immunotoxin,anti-(ricin A chain) immunotoxin

