



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## **Racotumomab Biosimilar - Anti-Idiotope of anti-(N-glycolylneuraminic acid (NeuGc,NGNA) -gangliosides GM3) Mus musculus IgM-kappa monoclonal antibody P3 mAb - Research Grade**

### **OVERVIEW**

---

<b>PRODUCT REFERENCE</b>	PX-TA1160
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1-kappa
<b>SOURCE</b>	CAS 946832-34-4

### **PRODUCT INFORMATION**

---

<b>ORIGIN SPECIES</b>	Mus musculus
<b>EXPRESSION SYSTEM</b>	Mammalian cells
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	Racotumomab,,Idiotope of anti-(N-glycolylneuraminic acid (NeuGc,NGNA) -gangliosides GM3) Mus musculus IgM-kappa monoclonal antibody P3,anti-Idiotope of anti-(N-glycolylneuraminic acid (NeuGc,NGNA) -gangliosides GM3) Mus musculus IgM-kappa monoclonal antib



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**