



Tremelimumab Biosimilar - Anti-CTLA4, CD152 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1177
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2-kappa
SOURCE	CAS 745013-59-6

PRODUCT INFORMATION

ORIGIN SPECIES	Homo sapiens
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Tremelimumab, Ticilimumab, CP-675, CP-675,206, CP-675206 clone 11.2.1, CTLA4, CD152, anti-CTLA4, CD152

