



Volociximab Biosimilar - Anti-ITGA5, CD49e mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1183
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-nd
SOURCE	CAS 558480-40-3

PRODUCT INFORMATION

ORIGIN SPECIES	Chimeric
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Volociximab,,ITGA5, CD49e,anti-ITGA5, CD49e

