



Tefibazumab Biosimilar - Anti-Fibrin-binding surface epitope clumping factor A mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1187
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
SOURCE	CAS 521079-87-8

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Tefibazumab,Aurexis,Fibrin-binding surface epitope clumping factor A,anti-Fibrin-binding surface epitope clumping factor A

