



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## Belimumab Biosimilar - Anti-TNFSF13B, CD257 mAb - Research Grade

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-TA1191
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1-lambda
<b>SOURCE</b>	CAS 356547-88-1

### PRODUCT INFORMATION

---

<b>ORIGIN SPECIES</b>	Homo sapiens
<b>EXPRESSION SYSTEM</b>	Mammalian cells
<b>BUFFER</b>	PBS buffer PH7.5
<b>STORAGE CONDITIONS</b>	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
<b>ALIAS/SYNONYMS</b>	Belimumab,LymphoStat-B,hBlySmAb-1.1,TNFSF13B, CD257,anti-TNFSF13B, CD257



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**