



Rovelizumab Biosimilar - Anti-ITGAL, CD11a, ITGB2, CD18 mAb - Research Grade

OVERVIEW

PRODUCT REFERENCE	PX-TA1216
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-nd
SOURCE	CAS 339086-79-2

PRODUCT INFORMATION

ORIGIN SPECIES	Humanized
EXPRESSION SYSTEM	Mammalian cells
BUFFER	PBS buffer PH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40 % glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.
ALIAS/SYNONYMS	Rovelizumab,Hu23F2G,ITGAL, CD11a, ITGB2, CD18,anti-ITGAL, CD11a, ITGB2, CD18

